
CHIROGUN

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DUO TRIGGERPUNKT DRÜCKER

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Stand: 10/2026

YOUR PRODUCT

The ChiroGun Duo Trigger Point Massager has two pressure points ideally spaced for the muscles to the left and right of the spine – the spine itself stays free in between. The suction cups hold it on smooth surfaces such as tiles or glass, so you can also use it standing against a wall.

USES

- For the back, left and right of the spine
- For shoulders, neck and glutes
- On the floor or on a wall via suction cups

TECHNICAL DATA

Size	approx. 14.5 cm long, balls approx. 5 cm in diameter each
Weight	approx. 230 g
Material	silicone
Mounting	2 suction cups
Colours	 Black  Red  Green

IN THE BOX

- 1 × ChiroGun Duo Trigger Point Massager in your chosen colour
- 1 × Storage pouch
- 1 × Exercise guide

HOW TO USE

1 Position

Place the Duo massager on a level floor or press the suction cups firmly onto a smooth, clean surface such as tiles, glass or a mirror. Check that it holds securely.

2 Hold pressure

Lie or lean against it so that your spine lies between the two balls. Hold the pressure for 20–60 seconds and keep breathing calmly.

3 Release

Slowly release the pressure. To detach it, gently lift the edge of a suction cup – do not yank the massager off.

KEY POINTS

- No more than 60 seconds per point, about 2–3 minutes per area.
- Always on the muscle, never directly on bones, joints or the spine.
- Control the pressure by supporting yourself with arms and legs.
- The pressure may be clearly felt, but never sharp or burning.

IMPORTANT

Read the safety information before first use. If you have any symptoms or pre-existing conditions, or are pregnant, talk to your doctor or physiotherapist first.

EXERCISES

EN


Always start with light pressure and increase it gradually.

1 UPPER BACK

Lie on your back, knees bent. Place the Duo massager crosswise under your upper back so that your spine lies between the two balls, and slowly let your weight sink onto it.

Keep your hips down for gentler pressure.

2 BACK AGAINST A WALL

Press the Duo massager's suction cups firmly onto a smooth wall at shoulder-blade height. Lean back against it so that your spine lies between the two balls, and slowly move up and down by bending your knees.

Check first that the suction cups hold securely.

CARE & STORAGE

- Wash with warm water and mild soap and leave to air-dry.
- Do not use solvents or harsh cleaning agents.
- Store protected from heat and direct sunlight.

DISPOSAL

Silicone does not go in the recycling collection. Please dispose of it in general waste unless your local authority specifies otherwise. Please separate the packaging by material for disposal.

WARRANTY

The statutory warranty applies. Please keep your proof of purchase. If you have any questions about the product, contact us at info@mineshop.de.

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mineshop.ch GmbH
mineShop GmbH