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DUOBALL 8 CM

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Stand: 10/2026

YOUR PRODUCT

The ChiroGun Duoball 8 cm consists of two connected balls. The groove between them keeps pressure off the spine: the balls massage the muscle strands to the left and right of it. With a diameter of 8 cm, it is ideal for the neck, upper back and soles of the feet.

USES

- For neck, upper back and shoulders
- For foot soles
- Small and light – ideal on the go

TECHNICAL DATA

Size	approx. 8 × 15 cm (diameter × length)
Weight	approx. 40 g
Material	100% EPP (expanded polypropylene), closed-cell structure
Properties	shape-stable, water-repellent, easy to clean
Temperature	shape-stable from -20 °C to +110 °C
Colours	 Black  Red  Green  Orange  Pink  Yellow

IN THE BOX

- 1 × ChiroGun Duoball 8 cm in your chosen colour

HOW TO USE

1 Prepare

Lay out a mat or find a level, non-slip spot. For wall exercises, you need a stable footing.

2 Massage

Place the Duoball across your spine so that your spine lies in the groove between the two balls. Slowly let your weight sink onto it and roll in small movements. Pause on tight spots for 10–30 seconds and keep breathing calmly.

3 Afterwards

Slowly ease off the Duoball and move around gently. Wipe it with a damp cloth if needed.

KEY POINTS

- Your spine always lies in the groove, never on a ball.
- About 1–2 minutes per area, no more than 30 seconds on one spot.
- Do not roll into your lower back and do not press on bones or joints.
- The pressure may be intense, but never sharp or burning.

IMPORTANT

Read the safety information before first use. If you have any symptoms or pre-existing conditions, or are pregnant, talk to your doctor or physiotherapist first.

EXERCISES

EN


The red arrows show the direction of movement. Always start with light pressure.

1 NECK

Lie on your back, knees bent. Place the Duoball crosswise under your neck with your neck vertebrae in the groove. Slowly turn your head left and right.


Gentle pressure only. Stop immediately if dizziness or tingling occurs.

2 UPPER BACK

Lie on your back, feet flat on the floor. Place the Duoball crosswise under your upper back with your spine in the groove. Cross your arms over your chest and roll slowly by lifting your hips.

Only roll to the middle of your back, not into your lower back.

3 BACK AGAINST A WALL

Stand with your back to a wall and place the Duoball crosswise between the wall and your back, with your spine in the groove. Roll slowly up and down by bending and straightening your knees.

How far your feet are from the wall sets the pressure.

4 SOLE OF THE FOOT

Standing, place one foot lengthwise on the Duoball and slowly roll from your heel to the ball of your foot.

Use a wall or chair back for support.

CARE & STORAGE

- Wipe with water and a little mild soap and leave to air-dry.
- Do not use harsh or abrasive cleaning agents or solvents.
- Store in a dry place, out of direct sunlight and not directly next to a heater, oven or open flame.
- Do not place heavy objects on the Duoball so that it keeps its shape.

DISPOSAL

The Duoball is made of polypropylene (PP, recycling code 05) and is recyclable. Dispose of it via your recycling collection or at a recycling centre and follow your local rules. Please separate the packaging by material for disposal.

WARRANTY

The statutory warranty applies. Please keep your proof of purchase. If you have any questions about the product, contact us at info@mineshop.de.

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mineshop.ch GmbH
mineShop GmbH