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
ERSATZ- BATTERIEN

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Stand: 10/2026

YOUR PRODUCT

The ChiroGun Spare Batteries are lithium-ion batteries for the ChiroGun Recovery Boots – one battery per boot. With a second set of charged batteries, the boots are always ready to use. The batteries are charged via USB-C, also separately from the boots.

USES

- As spare or additional batteries for the ChiroGun Recovery Boots
- For longer or several sessions a day

TECHNICAL DATA

Type	lithium-ion, 7.4 V per battery
Capacity	2,600 mAh
Charging	USB-C, approx. 2–3 hours
Size	approx. 8.5 cm long, approx. 4.5 cm diameter
Weight	approx. 250 g

IN THE BOX

- 1 × ChiroGun Battery Pack (2 batteries, one per boot)
- 1 × USB-C charging cable
- 1 × Instructions for safe use

HOW TO USE

1 Charge

Fully charge the batteries before first use – with the supplied USB-C charging cable. Charging time approx. 2–3 hours.

2 Insert

Switch off the boots, remove the empty batteries and insert the charged batteries into the boots' control units.

3 Store

Store batteries you are not using partially charged, dry and at room temperature.

KEY POINTS

- Only use for the ChiroGun Recovery Boots.
- Only charge with the supplied USB-C charging cable.
- If not used for a long time, recharge about every three months.
- Do not store fully discharged.

IMPORTANT

Read the safety information before first use. Lithium-ion batteries can overheat or catch fire if handled incorrectly.

CARE & STORAGE

- Only clean with a dry, soft cloth.
- Do not use water or cleaning agents.
- Store dry and protect from direct sunlight.

DISPOSAL

Batteries must not go in household waste (crossed-out wheeled bin symbol). Take used batteries to a battery collection point or a retailer. Tape over the contacts first to prevent short circuits.

WARRANTY

The statutory warranty applies. Please keep your proof of purchase. If you have any questions about the product, contact us at info@mineshop.de.

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mineshop.ch GmbH
mineShop GmbH