CHIROGUN

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Stand: 10/2026

YOUR PRODUCT

The ChiroGun Mini Pro is a compact cordless massage gun with a 24-watt motor and 4 intensity levels from gentle to strong. Weighing only 0.6 kg, it fits in any bag – ideal for on the go, the office or the gym.

USES

- To activate your muscles before training
- To loosen up after training or long days
- Compact for on the go, the office or travelling

TECHNICAL DATA

Motor	24 W
Settings	4 settings, 1,200–3,300 percussions/min
Noise	under 40 dB
Weight	approx. 0.6 kg
Charging	USB, approx. 2–3 hours
Battery life	up to 6 hours
Colours	 Black  Red  Green

IN THE BOX

- 1 × ChiroGun Mini Pro massage gun
- Massage attachments for different muscle groups
- 1 × USB charging cable
- 1 × Zip-up carry case

HOW TO USE

1 Charge

Fully charge the device before first use – only with the supplied charging cable or power adapter. Charge it on a firm, non-flammable surface and do not leave it unattended.

2 Choose an attachment

With the device switched off, push the right attachment firmly into place. The ball attachment suits almost all large muscle groups and is ideal to start with.

3 Massage

Switch the device on, start at the lowest setting and slowly move the massage head over the muscle. Do not press hard – the device does the work. 30–60 seconds per muscle group.

KEY POINTS

- Always start at the lowest setting and increase gradually.
- No more than 1–2 minutes on the same spot, about 15 minutes in total per session.
- Avoid bones, joints, spine, head, neck and face.
- Only change attachments when the device is off.

IMPORTANT

Read the safety information before first use. If you have a pacemaker or other implants or pre-existing conditions, or are pregnant, only use the device after consulting your doctor.

ATTACHMENTS


- Large, round attachments (e.g. ball) for large muscle groups such as thighs, calves and glutes.
- Flat attachments for an even massage of larger areas.
- Small, pointed or hard attachments only for individual points and with light pressure.

BATTERY AND CHARGING

- Only charge with the supplied charging cable or power adapter.
- Do not use the device while it is charging.
- Place on a firm, non-flammable surface for charging – not on a bed, sofa or under blankets.
- If not used for a long time, charge about every three months.
- Do not remove or open the battery yourself.

EXERCISES


Always start at the lowest setting and with light pressure.

1 FRONT OF THIGHS

Sit on a chair and, at a low setting, slowly move the massage head over your thigh between hip and knee.

Skip the kneecap. About 60 seconds per leg.

2 FOREARM

Rest your arm relaxed on the table and move the massage head between elbow and wrist with light pressure.

Avoid the wrist and elbow themselves.

CARE & STORAGE

- Switch off and unplug the charging cable before cleaning.
- Wipe the housing and attachments with a soft, slightly damp cloth – do not apply water directly to the device.
- Do not use corrosive cleaning agents or solvents.
- Store dry, at room temperature and in the supplied case.

DISPOSAL

Electrical devices and batteries must not go in household waste (crossed-out wheeled bin symbol). Take the device with its built-in battery to a collection point for waste electrical equipment or to a retailer that takes back old devices. Please separate the packaging by material for disposal.

WARRANTY

The statutory warranty applies. Please keep your proof of purchase. If you have any questions about the product, contact us at info@mineshop.de.

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mineshop.ch GmbH
mineShop GmbH