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
MINI ROLLE 15 CM

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Stand: 10/2026

YOUR PRODUCT

The ChiroGun Mini Roller is a small foam roller for targeted self-massage. With a length of 15 cm and a diameter of 6 cm, it reaches smaller areas such as the soles of the feet, forearms and neck. It is light and fits in any bag – ideal on the go and at work.

USES

- For feet, forearms, hands and neck
- As a short break at your desk
- On the go, at the gym or when travelling

TECHNICAL DATA

Size	approx. 15 × 6 cm (length × diameter)
Weight	approx. 20 g
Material	100% EPP (expanded polypropylene), closed-cell structure
Properties	shape-stable, water-repellent, easy to clean
Temperature	shape-stable from -20 °C to +110 °C
Colours	 Black  Red  Green  Orange  Pink  Yellow

IN THE BOX

- 1 × ChiroGun Mini Roller 15 cm in your chosen colour

HOW TO USE

1 Prepare

Place the mini roller on a level, non-slip surface – on the floor, a table or against a wall.

2 Roll

Place the mini roller under the desired spot and roll slowly back and forth with light pressure. Pause on tight spots for 10–20 seconds and keep breathing calmly.

3 Afterwards

Gently move the massaged area. Wipe the mini roller with a damp cloth if needed.

KEY POINTS

- About 1–2 minutes per area.
- Ideal for small areas such as feet, forearms and neck.
- Do not roll directly over bones, joints or the spine.
- The pressure may be intense, but never sharp or burning.

IMPORTANT

Read the safety information before first use. If you have any symptoms or pre-existing conditions, or are pregnant, talk to your doctor or physiotherapist first.

EXERCISES

EN


The red arrows show the direction of movement. Always start with light pressure.

1 SOLE OF THE FOOT

Sitting or standing, place one foot on the mini roller and slowly roll from your heel to the ball of your foot.

When standing, hold on to a wall or chair back.

2 FOREARM

Sit at a table, place the mini roller under your forearm and roll with light pressure between wrist and elbow.

Good as a short break when working at a computer.

3 NECK

Lie on your back, knees bent. Place the mini roller crosswise under your neck and slowly turn your head left and right.

Gentle pressure only. Stop immediately if dizzy.

CARE & STORAGE

- Wipe with water and a little mild soap and leave to air-dry.
- Do not use harsh or abrasive cleaning agents or solvents.
- Store in a dry place, out of direct sunlight and not directly next to a heater, oven or open flame.
- Do not place heavy objects on the mini roller so that it keeps its shape.

DISPOSAL

The roller is made of pure polypropylene (PP, recycling code 05) and is fully recyclable. Dispose of it via your recycling collection or at a recycling centre and follow your local rules. Please separate the packaging by material for disposal.

WARRANTY

The statutory warranty applies. Please keep your proof of purchase. If you have any questions about the product, contact us at info@mineshop.de.

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mineshop.ch GmbH
mineShop GmbH