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
NADELROLLER

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Stand: 10/2026

YOUR PRODUCT

The ChiroGun Needle Roller has a roller head with metal needles with slightly rounded tips and a wooden handle. It gently stimulates the skin and superficial fascia – to warm up before training or before a massage with a roller or ball.

USES

- To warm up before training
- To prepare for the foam roller or ball
- For arms, legs and back

TECHNICAL DATA

Size	approx. 26 × 9 × 6.5 cm
Weight	approx. 180 g
Needles	metal, slightly rounded tips
Handle	wood

IN THE BOX

- 1 × ChiroGun Needle Roller

HOW TO USE

1 Prepare

Your skin should be clean, dry and intact. Check the needle roller for bent or damaged needles before each use.

2 Roll

Roll the needle roller slowly and evenly over your arms, legs or back. Apply only light pressure – you should clearly feel the stimulation without irritating the skin. 2–5 minutes per area.

3 Afterwards

Clean the needle roller after each use and let it dry completely.

KEY POINTS

- Light pressure only – the needles must not injure the skin.
- 2–5 minutes per area, not several times a day on the same spot.
- Do not roll on the face or neck, or over bones or joints.
- For your personal use only – do not share, for hygiene reasons.


IMPORTANT

Read the safety information before first use. If you have any symptoms or pre-existing conditions, or are pregnant, talk to your doctor or physiotherapist first.

EXERCISES

EN


The red arrows show the direction of movement. Use light pressure only.

1 THIGHS

Sit on a chair and slowly roll the needle roller over your thigh between hip and knee with light pressure.

Skip the kneecap.

2 FOREARM

Rest your arm relaxed on the table and gently roll the needle roller between elbow and wrist with your other hand.

Skip the wrist and elbow.

CARE & STORAGE

- After each use, clean with warm water, soap and a soft brush.
- Dry thoroughly, e.g. in the air or with a hairdryer on the cold setting.
- Do not leave the wooden handle in water for long.
- Store dry and out of reach of children.

DISPOSAL

Please dispose of it in general waste unless your local authority specifies otherwise. Wrap the needles so that nobody can get hurt. Please separate the packaging by material for disposal.

WARRANTY

The statutory warranty applies. Please keep your proof of purchase. If you have any questions about the product, contact us at info@mineshop.de.

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mineshop.ch GmbH
mineShop GmbH