CHIROGUN NADELROLLER

Sicherheitsinformationen · Safety information · Consignes de sécurité · Informazioni di sicurezza ·
Información de seguridad · Veiligheidsinformatie · Informacje o bezpieczeństwie · Bezpečnostní
informace · Bezpečnostné informácie

DE	Deutsch	Sicherheitsinformationen	2
EN	English	Safety information	4
FR	Français	Consignes de sécurité	6
IT	Italiano	Informazioni di sicurezza	8
ES	Español	Información de seguridad	10
NL	Nederlands	Veiligheidsinformatie	12
PL	Polski	Informacje o bezpieczeństwie	14
CS	Čeština	Bezpečnostní informace	16
SK	Slovenčina	Bezpečnostné informácie	18

Variante / Variant	Artikelnummer / Item number	EAN
–	10203124	7649995881046

SICHERHEITSINFORMATIONEN

ChiroGun Nadelroller – Nadelroller zur Hautstimulation

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SAFETY INFORMATION

ChiroGun Needle Roller – Needle roller for skin stimulation

Please read completely before first use and keep with the product. Pass this information on if you give the product to someone else.

Intended use

- For gentle skin stimulation on the arms, legs and back of adults.
- Not a medical device. Does not replace treatment by a doctor or physiotherapist.
- For personal use only – do not share with others, for hygiene reasons.

Do not use – or only after consulting a doctor – if you

- have open wounds, injuries, inflammation, sunburn or irritated skin
- have skin conditions such as eczema, psoriasis, acne or infections
- have moles, warts or scars at the site
- have varicose veins, vein disorders or a blood clotting disorder, or take blood-thinning medication
- have diabetes with impaired sensation or other pre-existing conditions, or are pregnant

During use

- Do not roll on the face, neck, mucous membranes or intimate area, or over bones and joints.
- Apply light pressure only. Stop immediately in case of redness, burning, pain or slight bleeding.
- Bent or damaged needles: stop using the product.

Children

- Not a toy. Keep out of reach of children – risk of injury from sharp needles.

Hygiene

- Clean before and after each use and let it dry completely.

MATERIAL

Metal needles with slightly rounded tips, wooden handle.

MANUFACTURER AND RESPONSIBLE PERSON

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This safety information belongs with the product. Please keep it. Version: 10/2026

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mineshop.ch GmbH
mineShop GmbH