# CHIROGUN

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# REGENERATIONS STIEFEL

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Stand: 10/2026

# YOUR PRODUCT

The ChiroGun Recovery Boots are cordless compression boots for the legs. Air chambers build up pressure in cycles of about 60 seconds and then release it. The compressors sit directly on the boot and run on removable batteries – ideal for recovery after sport, before training or after long days of sitting or standing.

# USES

- For leg recovery after sport
- Before training or after sitting or standing for long
- At home on the sofa or on the go

# TECHNICAL DATA

|                   |                                                                        |
|-------------------|------------------------------------------------------------------------|
| Sizes             | S: approx. 84 × 33 cm · M: approx. 94 × 38 cm · L: approx. 102 × 38 cm |
| Weight            | per boot incl. battery: S approx. 1,160 g, M/L approx. 1,370 g         |
| Pressure          | 20–120 mmHg                                                            |
| Battery           | lithium-ion, removable, 1,000 mAh, 7.4 V, 7.4 Wh                       |
| Charging          | USB-C, approx. 1.5 hours, 5 V/2 A power adapter (not included)         |
| Run time          | approx. 60 minutes at medium pressure                                  |
| Power / noise     | 25 W / up to 55 dB                                                     |
| Protection rating | IP20                                                                   |
| Material          | polyamide with TPU coating, YKK zips                                   |

# IN THE BOX

- 2 × ChiroGun Recovery Boot
- 2 × Battery in a battery pouch
- 1 × USB-C dual charging cable
- 1 × Storage bag made of recycled PET
- 1 × Quick guide

# HOW TO USE

## 1 Charge

Fully charge both batteries before first use – with the supplied USB-C charging cable and a 5 V/2 A power adapter (not included). Charging time approx. 1.5 hours.

## 2 Put on

Sit or lie down comfortably with your feet no lower than your pelvis. Put on the boots, close the zips and insert the charged batteries with the control units.

## 3 Use

Switch the boots on and start with low pressure. After exercise, 30–60 minutes at medium to high intensity; before exercise, 15–30 minutes at low to medium intensity.

### KEY POINTS

- Always start with low pressure and increase gradually.
- Only use lying down or sitting, not standing or walking.
- The pressure may feel firm, but never painful. Switch off immediately if you feel tingling or numbness.
- No more than one 60-minute session in one go.

### IMPORTANT

Read the safety information before first use. If you have vein, heart or vascular disease, diabetes or a pacemaker, or are pregnant, only use the boots after consulting your doctor.

## CHOOSE YOUR SIZE


- Thigh circumference under 67 cm: size S; under 75 cm: size M or L.
- Inside leg under 79 cm: size S; 80–86 cm: size M.
- How to measure your inside leg: stand upright against a wall, clamp a book between your legs at the crotch and measure from the floor to the top edge of the book.

## BATTERY AND CHARGING

- Only charge with the supplied USB-C charging cable and a 5 V/2 A power adapter.
- Do not use the boots while they are charging.
- Place on a firm, non-flammable surface to charge and do not leave unattended.
- The batteries are removable and can be charged separately.

## POSITION FOR USE

EN

### 1 CORRECT POSITION

Sit or lie comfortably, legs stretched out, feet no lower than your pelvis. This way the compression works with gravity.

*Ideal on a sofa or bed with a backrest.*

## CARE & STORAGE

- Switch off and remove the batteries before cleaning.
- Wipe with a slightly damp cloth – do not apply water directly to the control unit or batteries.
- Do not machine wash. Do not use corrosive cleaning agents or solvents.
- Store dry in the storage bag; do not kink or place heavy objects on them.

## DISPOSAL

Electrical devices and batteries must not go in household waste (crossed-out wheeled bin symbol). Take the boots, control units and batteries to a collection point for waste electrical equipment or to a retailer that takes back old devices. Please separate the packaging by material for disposal.

## WARRANTY

The statutory warranty applies. Please keep your proof of purchase. If you have any questions about the product, contact us at [info@mineshop.de](mailto:info@mineshop.de).

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mineshop.ch GmbH
mineShop GmbH