CHIROGUN REGENERATIONS STIEFEL

Sicherheitsinformationen · Safety information · Consignes de sécurité · Informazioni di sicurezza ·
Información de seguridad · Veiligheidsinformatie · Informacje o bezpieczeństwie · Bezpečnostní
informace · Bezpečnostné informácie

DE	Deutsch	Sicherheitsinformationen	2
EN	English	Safety information	5
FR	Français	Consignes de sécurité	8
IT	Italiano	Informazioni di sicurezza	11
ES	Español	Información de seguridad	14
NL	Nederlands	Veiligheidsinformatie	17
PL	Polski	Informacje o bezpieczeństwie	20
CS	Čeština	Bezpečnostní informace	23
SK	Slovenčina	Bezpečnostné informácie	26

Größe / Size	Artikelnummer / Item number	EAN
S	10203091	7649995881305
M	10203092	7649995881299
L	10203093	7649995881282

SICHERHEITSINFORMATIONEN

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SAFETY INFORMATION

ChiroGun Recovery Boots – Compression boots for the legs

Please read completely before first use and keep with the product. Pass this information on if you give the product to someone else.

Intended use

- For compression massage of adults' legs, for private use, while sitting or lying down.
- Not a medical device. Does not replace treatment by a doctor or physiotherapist, nor lymphatic drainage.
- Not to be used by children. Persons with reduced physical, sensory or mental capabilities only under supervision.

Do not use – or only after consulting a doctor – if you

- have (suspected) thrombosis, vein inflammation or varicose veins
- have heart failure, cardiovascular disease, high blood pressure or fluid in the lungs
- have circulatory disorders in the legs (e.g. PAD) or diabetes with impaired sensation
- have open wounds, skin inflammation, infections, recent injuries or fractures on your legs
- have a pacemaker or other electronic implant
- are pregnant, have tumours or have had surgery

During use

- Only use sitting or lying down – not while standing, walking or sleeping.
- Start with low pressure. Switch off immediately in case of pain, tingling, numbness, skin discolouration or dizziness.
- Use for no more than 60 minutes at a time.
- Do not put on over jewellery, zips or hard objects on clothing.

Electrical safety and battery

- Only charge with the supplied charging cable and a 5 V/2 A power adapter. Do not use while charging.
- Keep the control unit and batteries away from water. Do not operate with wet hands; do not use in the bathroom.
- Do not use if the boots, hoses, cables, control unit or batteries are damaged, or if they get hot, smoke or smell unusual.
- Do not open or repair them yourself. Do not short-circuit or puncture the batteries or expose them to high heat.

Children

- Not a toy. Keep the boots, batteries and cables out of reach of children – risk of strangulation from cables.

Heat and fire

- Do not charge or store near heat sources or open flames, or in direct sunlight.

MATERIAL

Polyamide (PA) with TPU coating, YKK zips, bag made of 100% PET. Control units with removable lithium-ion batteries.

MANUFACTURER AND RESPONSIBLE PERSON

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This safety information belongs with the product. Please keep it. Version: 10/2026

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mineshop.ch GmbH
mineShop GmbH