
CHIROGUN

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Stand: 10/2026

YOUR PRODUCT

The ChiroGun Trigger Point Massager is a massage tool for targeted, pinpoint pressure. You lie or lean against the tip and let your body weight slowly work on tense spots – on the floor or against a wall, e.g. on the shoulders, back and glutes.

USES

- For shoulders, upper back and glutes
- On the floor, against a wall or seated
- Small and light – ideal on the go

TECHNICAL DATA

Size	approx. 7.5 cm high, approx. 7 cm diameter
Weight	approx. 150 g
Material	silicone, free of latex, BPA and plasticisers (per manufacturer)
Colours	 Black  Green  Red

IN THE BOX

- 1 × ChiroGun Trigger Point Massager in your chosen colour
- 1 × Quick guide

HOW TO USE

1 Position

Place the trigger point massager on a level, firm surface or hold it against a wall. Lie or lean against it so that the tip meets the muscle you want to treat.

2 Hold pressure

Slowly let your body weight sink onto the point and hold the pressure for 20–60 seconds. Keep breathing calmly and deeply.

3 Release

Slowly release the pressure, gently move the area and find the next point if needed.

KEY POINTS

- No more than 60 seconds per point, about 2–3 minutes per area.
- Always on the muscle, never directly on bones, joints or the spine.
- Control the pressure by supporting yourself with arms and legs.
- The pressure may be clearly felt, but never sharp or burning.

IMPORTANT

Read the safety information before first use. If you have any symptoms or pre-existing conditions, or are pregnant, talk to your doctor or physiotherapist first.

EXERCISES

EN

Always start with light pressure and increase it gradually.

1 UPPER BACK


Lie on your back, knees bent. Place the massager next to your spine, under the muscle between your shoulder blades, and slowly let your weight sink onto it.

Keep your hips down for gentler pressure.

2 GLUTES

Sitting, support yourself with your hands behind you and place the massager under one buttock. Adjust the pressure by shifting your weight. Then switch sides.

3 SHOULDER AGAINST A WALL

Hold the massager between the wall and your shoulder blade and slowly lean against it. Move the point by bending your knees slightly.

Against a wall, the pressure is easiest to control.

CARE & STORAGE

- Wash with warm water and mild soap and leave to air-dry.
- Do not use solvents or harsh cleaning agents.
- Store protected from heat and direct sunlight.

DISPOSAL

Silicone does not go in the recycling collection. Please dispose of it in general waste unless your local authority specifies otherwise. Please separate the packaging by material for disposal.

WARRANTY

The statutory warranty applies. Please keep your proof of purchase. If you have any questions about the product, contact us at info@mineshop.de.

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mineshop.ch GmbH
mineShop GmbH